



## Technical Datasheet

<b>Product:</b>	10µl Short Pipette Tips
<b>Part No:</b>	310-xxxx
<b>Comments:</b>	Universal Fit



### Specification

<b>Product:</b>	<ul style="list-style-type: none"><li>• 10µl Short Pipette Tips</li></ul>
<b>Dimension, Length:</b>	<ul style="list-style-type: none"><li>• 29mm</li></ul>
<b>Colour, Optical Quality:</b>	<ul style="list-style-type: none"><li>• Natural, virgin, PP, translucent</li></ul>
<b>Low Retention:</b>	<ul style="list-style-type: none"><li>• Low adhesion, ultra-hydrophobic surface PP 310-4050</li></ul>
<b>Sterility:</b>	<ul style="list-style-type: none"><li>• By e-Beam, 310-4050; 310-4150</li></ul>
<b>Graduation:</b>	<ul style="list-style-type: none"><li>• At 2µl and 10µl</li></ul>
<b>Nominal Capacity:</b>	<ul style="list-style-type: none"><li>• 10µl</li></ul>
<b>Material / Characteristics:</b>	<ul style="list-style-type: none"><li>• Natural Polypropylene (PP), free of heavy metal and softener</li><li>• Autoclavable at maximum 121°C</li><li>• Use temperature -20°C up to 90°C</li></ul>
<b>Product Characteristics:</b>	<ul style="list-style-type: none"><li>• Physiological harmless / non-toxic</li><li>• Biologically inert</li><li>• RNase, DNase, DNA and Pyrogen-free</li><li>• PrecisePointe, features thin walls and small orifice</li></ul>
<b>Packing, Bulk:</b>	<ul style="list-style-type: none"><li>• 1000 Tips/Bag, 10 Bags/Case</li><li>• Case specifications: 26 x 26 x 18cm; 0.01cu.m; 1.36kg</li></ul>
<b>Packing, Reload:</b>	<ul style="list-style-type: none"><li>• 10x 96 Tips/Rack, 10x Racks/Pack, 10x Packs/Case</li><li>• Case specifications: 39 x 31 x 30cm; 0.04cu.m; 5.44kg</li></ul>
<b>Packing, Racked:</b>	<ul style="list-style-type: none"><li>• 10x 96 Tips/Rack, 10x Racks/Pack. 5x Packs/Case</li><li>• Case specifications: 39 x 31 x 30cm; 0.04cu.m; 5.90kg</li></ul>
<b>Storage:</b>	<ul style="list-style-type: none"><li>• In dry condition</li><li>• Double boxes and master shipper cases originally closed</li><li>• Storage temperature +4°C to maximal 25 °C</li><li>• Pallets should not be stacked</li></ul>
<b>Conformities / Regulations:</b>	<ul style="list-style-type: none"><li>• Rack conforms to the SBS footprint/configuration compatible with robotic workstations, multichannel pipettors, and auto-sampling instruments</li><li>• Compatibility - fits Gilson Pipetman P2/P10 and Eppendorf Ultra Micro style. For full compatibility see BioPointe table.</li></ul>