



## Technical Datasheet

<b>Product:</b>	MicroTube Rack System (sometimes known as Cluster Tubes)
<b>Part No:</b>	STN3346-12C
<b>Comments:</b>	Deep-well microtiter plates with removable wells, available sterile

### Specification

<b>Product:</b>	<ul style="list-style-type: none"><li>• MicroTube Rack System, 12-CapStrips, Vented in bulk packs</li></ul>
<b>Colour, Optical Quality:</b>	<ul style="list-style-type: none"><li>• Natural, virgin, PE, translucent</li></ul>
<b>Graduation:</b>	<ul style="list-style-type: none"><li>• n/a</li></ul>
<b>Nominal Capacity:</b>	<ul style="list-style-type: none"><li>• n/a</li></ul>
<b>Material / Characteristics:</b>	<ul style="list-style-type: none"><li>• Natural Polyethylene (PE), free of heavy metal or additives</li><li>• Not autoclavable</li><li>• Use temperature -20°C up to 90°C</li></ul>
<b>Product Characteristics:</b>	<ul style="list-style-type: none"><li>• Biologically inert</li><li>• Physiologically harmless / non-toxic</li><li>• RNase, DNaseI, DNA and Pyrogen-free</li></ul>
<b>Packing:</b>	<ul style="list-style-type: none"><li>• 80 CapStrips/Bag</li><li>• 5 Bags/Case</li><li>• Case specifications: 20.4 x 15.4 x 11.9cm; 0.91kg</li></ul>
<b>Storage:</b>	<ul style="list-style-type: none"><li>• In dry condition</li><li>• Storage temperature +4°C to maximal 25 °C</li><li>• Pallets should not be stacked</li></ul>
<b>Conformities / Regulations:</b>	<ul style="list-style-type: none"><li>• Conforms to the SBS footprint/configuration compatible with robotic workstations, multichannel pipettors, and auto-sampling instruments</li></ul>

